

# T-Edge Autoclave

The only gravity autoclave  
upgradable to pre/post vacuum  
by software!

- Larger chamber capacity and shorter cycle times
- Convenient and easy loading with door opening to 110°
- Available in 110V and 220V
- 10" and 11" Chambers



Designed to meet all current  
sterilization standards ANSI/AAMI ST55

  
**Tuttnauer**  
Innovation · Legacy · Partnership

## Features:

- Optimal chamber capacity
- Convenient and easy loading with door opening to 110°
- Assorted water filling and drainage methods (top, front or automatic)
- Range of automatic one-button programs
- User-friendly touchpad and control panel
- Advanced documentation and traceability options
- Remote Wi-Fi monitoring via smartphone or tablet

Printer includes sterilization data (including expiration date).



## Technical Specifications

Autoclave Model	Chamber Size (Inches)	Volume	Overall Dimensions D x W x H (Inches)	Number of trays	Connections	Standard Cassette Capacity	Item Supplies	Power Supply	Power Consumption	Weight
T-Edge 10 220V	18" depth x 10" diameter	6 gal/ 23 liters	21.1"x18.89" x19.68"	5	1 ethernet port, 4 USB ports, WiFi, intergrated automatic water filling and drain connection	3 Full + 3 Half	Drain tube, tray holder	230V, Frq 50-60Hz	2000W	110 lbs.
T-Edge10 110V	18" depth x 10" diameter	6 gal/ 23 liters	21.1"x18.89" x19.68"	5	1 ethernet port, 4 USB ports, WiFi, intergrated automatic water filling and drain connection	3 Full + 3 Half	Drain tube, tray holder	120V, Frq 50-60Hz	1700W	110 lbs.
T-Edge11 110V	18" depth x 11" diameter	7 gal/ 27 liters	21.1"x19.68" x19.68"	5	1 ethernet port, 4 USB ports, WiFi, intergrated automatic water filling and drain connection	4 Full + 4 Half	Drain tube, tray holder	120V, Frq 50-60Hz	1700W	115 lbs.

## Sterilization Programs

Class S	Icon	Cycle Program	Load Type	*T-Edge10 Cycle Times Including Drying (Minutes)	*T-Edge11 Cycle Times Including Drying (Minutes)
		Unwrapped Instruments 270° F	Unwrapped Instruments	20 Minutes	25 Minutes
		Wrapped Pouches 270° F	Wrapped Pouches	32 Minutes	55 Minutes
		Unwrapped Delicate 250° F	Unwrapped Delicate Instruments	45 Minutes	47 Minutes
		Handpieces 270° F	Handpieces	45 Minutes	60 Minutes

Class B	Icon	Cycle Program	Load Type	*T-Edge10 Cycle Times Including Drying (Minutes)	*T-Edge11 Cycle Times Including Drying (Minutes)
		Unwrapped Instruments 273° F	Unwrapped Instruments	20 Minutes	26 Minutes
		Wrapped Pouches 273° F	Wrapped Pouches	39 Minutes	56 Minutes
		Unwrapped Delicate 250° F	Unwrapped Delicate Instruments	35 Minutes	35 Minutes
		Wrapped Delicate 250° F	Wrapped Delicate Instruments	60 Minutes	79 Minutes

### Sterisil® Ac+ System

Convenient, on demand autoclave water, made right in your practice



STERISIL23

## Tuttnauer

Innovation · Legacy · Partnership

Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720

E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

Learn from our experts.

Visit our blog:

<https://tuttnauerusa.com/blog>

